





PRODUCT SPECIFICATIONS

Date: 2013-09-20

Note: For research use only

PRODUCT NAME : Anti-HBsAg 2506 SPTN-5

PRODUCT SPECIFICITY : Antibody recognizes Hepatitis B virus surface antigen

PRODUCT CODE : 100353

PRODUCT BUFFER : 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.1 % NaN $_3$ as a preservative

SHELF LIFE AND STORAGE : Shelf life unspecified, storage at 2-8 °C

ANALYTE DESCRIPTION : Hepatitis B is an infectious illness caused by hepatitis B virus (HBV). The acute

illness causes liver inflammation, vomiting, jaundice and rarely, death. Chronic hepatitis B may eventually cause liver cirrhosis and liver cancer—a fatal disease with very poor response to current chemotherapy. The infection is preventable by vaccination. The hepatitis B surface antigen (HBsAg) is most frequently used to screen for the presence of this infection. It is the first detectable viral antigen to appear during infection. However, early in an infection, this antigen may not be present and it may be undetectable later in the infection as it is being cleared by the

nost

PARAMETERS TESTED FROM EACH LOT

PRODUCT APPEARANCE : Clear liquid, may turn opaque during storage

PRODUCT ACTIVITY : 80 – 120 % compared to reference in an IFMA-test

IEF-RANGE : 8.4 – 8.6

PURITY : ≥ 95 %

PARAMETERS DETERMINED ONLY DURING PRODUCT R&D PHASE

CLASS AND SUBCLASS : IgG2a

ASSOCIATION CONSTANT : N/D

DISSOCIATION CONSTANT : N/D

AFFINITY CONSTANT : N/D

DETERMINATION METHOD : N/A

ANTIGEN : N/A

CROSS-REACTIVITIES : Not tested



INFECTIOUS DISEASES

medixMar[™]

EPITOPE : N/D

EPITOPE GROUP : N/A

Two antibodies binding to the same, or closely located epitopes, belong to the same group and hence cannot be used as a pair in a sandwich assay. Epitope group numbering does not give any detailed information where the epitope is located.

PAIR RECOMMENDATIONS : SOLID LABEL

Please note that pair recommendations are based on results obtained in our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations should be taken only as a directive.

PRODUCT STABILITY

Please note that the shelf life given on page one is based on real time stability testing at +2-8 °C in the product buffer.

TEMPERATURE, DAYS	RESULT
-70 °C, 21 days	N/D
-20 °C, 21 days	N/D
+4 °C, 21 days	N/D
+35 °C, 7 days	N/D
+35 °C, 21 days	N/D
+45 °C, 3 days	N/D
+45 °C, 7 days	N/D

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. The maximum duration of the test is 21 days, except for the $+45\,^{\circ}\text{C}$ only 7 days.

pH, 14 DAYS, +4°C	RESULT
5.0	N/D
6.0	N/D
7.0	N/D
8.0	N/D

Stability testing is performed to see whether pH affects the antigen binding, charge or composition of the antibody during 14 days at +4 °C.

MISCELLANEOUS : -

REFERENCES :

Legal disclaimer

MedixMAB monoclonal antibodies meet their specifications if transported, stored and used according to the instructions for use. See also separate Notes on use and storage of monoclonal antibodies. If not otherwise agreed in writing all products are sold under General Delivery Terms. Other terms may apply if MedixMAB monoclonal antibodies are purchased from an official distributor of Medix Biochemica products. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission. Further information on products and methods is available from medixMAB website.