



PRODUCT SPECIFICATION

catalog number:	2-024-C
product description:	Cardiac Troponin T Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Highly purified human cardiac troponin T subunit
form:	Liquid
purification method:	Protein A chromatography
purity:	>95% by SDS-PAGE
buffer:	10 mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2
preservative:	0.05 % (w/v) NaN ₃
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.