



## PRODUCT SPECIFICATION

catalog number:	<b>A10751</b>
product description:	NT-proBNP (38-44) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Native human proBNP
source:	Murine
form:	Liquid
purity:	>90% pure
buffer:	10mM PBS, pH 7.4
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific for the 38-44aa sequence of NT-proBNP.
pairing:	This antibody can be used in ELISA as the capture antibody in conjunction with A10741 or A10742 (as the conjugate). This clone can also be mixed with monoclonal A10752.
expiration date:	[lot specific]
storage:	2-8°C short term (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.