



PRODUCT SPECIFICATION

catalog number:	A91064
product description:	Procalcitonin (PCT) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purity:	>90%
buffer:	PBS
preservative:	0.1% Proclin 300
isotype:	IgG
concentration:	[lot specific] mg/mL
pairing:	This antibody can be used as the detection antibody in combination with A91062 as the capture antibody.
expiration date:	[lot specific]
storage:	2-8°C short term (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.