



PRODUCT SPECIFICATION

catalog number:	D22510								
product description:	DHEA-Sulfate Polyclonal Antibody								
lot number:	[lot specific]								
immunogen:	DHEA-3-HS-BSA								
source:	Rabbit								
form:	Neat antiserum								
preservative:	None								
titer:	[lot specific] tests per mL (by RIA)								
cross-reactivity:	<table><tr><td>▪ DHEA-Sulfate</td><td>100%</td><td>▪ Other related steroids</td><td><1%</td></tr><tr><td>▪ DHEA</td><td>50%</td><td></td><td></td></tr></table>	▪ DHEA-Sulfate	100%	▪ Other related steroids	<1%	▪ DHEA	50%		
▪ DHEA-Sulfate	100%	▪ Other related steroids	<1%						
▪ DHEA	50%								
expiration date:	[lot specific]								
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.								
hazard/biohazard:	There is no known hazard or biohazard associated with this product.								