



PRODUCT SPECIFICATION

CATALOG NUMBER: D23010

PRODUCT DESCRIPTION:	Rabbit anti Progesterone -11 HSA
HOST ANIMAL:	Rabbit
IMMUNOGEN:	Progesterone-11-HSA
TITER:	12,000,000 RIA
STORAGE TEMPERATURE:	2-8°C short term -20° C long term
CROSS REACTIVITY:	Progesterone 11 100% 17-alpha-Hydroxyprogesterone <0.1% DeoxyCorticosterone <1.7% Corticosterone <2.0% 20-alpha-hydroxyprogesterone < 0.2%
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.