



PRODUCT SPECIFICATION

CATALOG NUMBER: G05010

PRODUCT DESCRIPTION:	Rabbit anti NGAL (Lipcalin-2), Affinity Purified
FORM:	Liquid
PURITY:	>95% by Caprylic Acid Ammonium Sulfate Precipitation
BUFFER:	0.02M HEPES, pH 8.0
SUBCLASS:	IgG
STORAGE:	-20°C or -80°C. Avoid repeated freeze/thaw cycles.
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this material.