



PRODUCT SPECIFICATION

CATALOG NUMBER: J05010

PRODUCT DESCRIPTION:	Recombinant NGAL (Mutant) Antigen
SOURCE:	<i>E.coli</i>
FORM:	Liquid
PURITY:	> 95% by SDS-PAGE
PURIFICATION:	Column Chromatography
BUFFER:	50mM sodium phosphate, 0.15M sodium chloride, 50% glycerol, pH 7.4 ± 0.2
BIOBURDEN:	< 10 CFU/ml
STORAGE:	-20°C; avoid repeat freeze-thaw cycles
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product