



PRODUCT SPECIFICATION

CATALOG NUMBER: J05011

PRODUCT DESCRIPTION:	Recombinant NGAL (Wild Type) Antigen
SOURCE:	<i>E.coli</i>
FORM:	Liquid
PURITY:	>95% by SDS-PAGE ^{12%}
PURIFICATION:	Column Chromatography
BUFFER:	50mM sodium phosphate, 0.15M sodium chloride, 50% glycerol, pH 7.4 ± 0.2
BIOBURDEN:	< 10 CFU/ml
STORAGE:	-20°C; avoid repeat freeze-thaw cycles
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product