



## PRODUCT SPECIFICATION

**CATALOG NUMBER: 1-023**

PRODUCT DESCRIPTION:	Purified Human Cardiac Myoglobin Antigen
SOURCE:	Human heart
FORM:	Liquid
PURITY:	>98% by SDS-PAGE
BUFFER:	50 mM Tris, pH 8.0
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE TEMPERATURE:	2-8°C
HAZARD/BIOHAZARD:	Donor's serum was tested and found negative for HbSAg, HIV-1 antigen, HIV-1, HIV-2, and HCV antibodies by FDA approved tests. This material contains sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.