



PRODUCT SPECIFICATION

catalog number:	J46833
product description:	Zika NS-1HEK Antigen, Recombinant, 6xHis tag at C-terminus
lot number:	[lot specific]
source:	<i>Human Embryonic Kidney 293 cells</i>
form:	Liquid, sterile filtered
purity:	>90% by SDS-PAGE
buffer:	Dulbecco's Phosphate buffer saline, pH7.4
concentration:	[lot specific] mg/mL
storage:	2-8°C short term (up to 1 week); -18°C or lower long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known biohazard associated with this product.