



PRODUCT SPECIFICATION

catalog number:	V24100
product description:	TCA-BSA Conjugate
lot number:	[lot specific]
form:	Liquid
buffer:	PBS, pH 7.2
preservative:	0.2% NaN ₃
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.