



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>V50042</b>   |
| product description: | K2-BSA Conjugate  |
| lot number:          | [lot specific]  |
| form:                | Liquid  |
| buffer:              | Distilled water   |
| preservative:        | 0.01% NaN <sub>3</sub>  |
| concentration:       | [lot specific] mg/mL  |
| hapten/protein:      | 30:1  |
| pairing:             | This conjugate can be paired with K2 monoclonal antibody A50042.  |
| notes/applications:  | This conjugate is known to work well in lateral flow.   |
| storage:             | 2-8°C (≤ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.  |
| hazard/biohazard:    | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |