



## PRODUCT SPECIFICATION

catalog number:	<b>V50050</b>
product description:	MDPV-BSA Conjugate
lot number:	[lot specific]
form:	Liquid
buffer:	PBS
preservative:	none
concentration:	[lot specific] mg/mL
hapten/protein:	17:1
pairing:	This conjugate can be paired with MDPV monoclonal antibody A50055.
notes/applications:	This conjugate is known to work well in lateral flow.
storage:	2-8°C short term ( $\leq$ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.