



## PRODUCT SPECIFICATION

catalog number:	<b>V55050</b>
product description:	Phencyclidine (PCP)-BSA Conjugate
lot number:	[lot specific]
form:	Liquid
buffer:	PBS, pH 7.2
preservative:	0.02% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
expiration date:	This conjugate is known to be stable for six months at 2-8°C or up to one year at -20°C from the date of receipt.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.