



PRODUCT SPECIFICATION

catalog number:	V56046
product description:	Ethyl Glucuronide (EtG)-BSA Conjugate
lot number:	[lot specific]
form:	Liquid
buffer:	Distilled water
preservative:	0.01% NaN ₃
concentration:	[lot specific] mg/mL
hapten/protein:	16:1
pairing:	This conjugate can be paired with EtG monoclonal antibody A56044.
notes/applications:	This conjugate is known to work well in lateral flow and ELISA.
storage:	2-8°C short term (≤ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.